



CDCA2 Antibody

Product Code	CSB-PA005025GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q69YH5
Immunogen	Human CDCA2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	cell division cycle associated 2;CDCA2;FLJ25804;MGC129906;MGC129907;Repo-Man ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDCA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.