



CENPL Antibody

Product Code	CSB-PA005215GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0S6
Immunogen	Human CENPL
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	centromere protein L;CENPL;C1orf155;CENP-L;FLJ31044;FLJ31786;RP3-383J4.1;dJ383J4.3 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CENPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.