



Recombinant Prosthecochloris vibrioformis Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033PZJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4SG02
Product Type	Recombinant Protein
Immunogen Species	Chlorobium phaeovibrioides (strain DSM 265 / 1930) (Prosthecochloris vibrioformis (strain DSM 265))
Purity	>85% (SDS-PAGE)
Sequence	MFTPYGMSTL VKLAAVSVGT LLFSILLPLQ AAIPVGSAAA IFLIFSLWFF RDPNRS LPPG KGIV IAPADG TVMSIKACSH PFTGNGSSII SIFMSPLNVH VNRIPVTGTV TLLQHHPGSF SMAFDEKSGE ENERMEIGIE SNGMKLHFTQ VAGFLARRIV CPLSLNEPVT AGKRFGMIKF G
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.