



CHCHD6 Antibody

Product Code	CSB-PA005321GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRQ6
Immunogen	Human CHCHD6
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	coiled-coil-helix-coiled-coil-helix domain containing 6; CHCHD6; MGC13016 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHCHD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.