



# CHD4 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA005329GA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | Q14839   |
| <b>Immunogen</b>           | Human CHD4   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human,Mouse,Rat  |
| <b>Tested Applications</b> | ELISA,WB,IF  |
| <b>Storage Buffer</b>      | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.  |
| <b>Purification Method</b> | Antigen Affinity purified  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | chromodomain helicase DNA binding protein<br>4;CHD4;DKFZp686E06161;Mi-2b;Mi2-BETA ;  |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody   |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | CHD4   |
| <b>Target Details</b>      | The product of this gene belongs to the SNF2/RAD54 helicase family. It represents the main component of the nucleosome remodeling and deacetylase complex and plays an important role in epigenetic transcriptional repression. Patients with dermatomyositis develop antibodies against this protein. |