



CHMP1B Antibody

Product Code	CSB-PA005359GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7LBR1
Immunogen	Human CHMP1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chromatin modifying protein 1B;CHMP1B;C10orf2;C18-ORF2;C18orf2;CHMP1.5;Vps46-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHMP1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.