



# CHST14 Antibody

|                            |                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA005407GA01HU                                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                   |
| <b>Uniprot No.</b>         | Q8NCH0                                                                                          |
| <b>Immunogen</b>           | Human CHST14                                                                                    |
| <b>Raised In</b>           | Rabbit                                                                                          |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                 |
| <b>Tested Applications</b> | ELISA,WB,IHC                                                                                    |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.            |
| <b>Purification Method</b> | Antigen Affinity purified                                                                       |
| <b>Isotype</b>             | IgG                                                                                             |
| <b>Alias</b>               | carbohydrate (N-acetylgalactosamine 4-0) sulfotransferase 14;CHST14;D4ST-1;D4ST1;HD4ST;HNK1ST ; |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                                  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                            |
| <b>Target Names</b>        | CHST14                                                                                          |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                     |