



CHTF18 Antibody

Product Code	CSB-PA005418GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVB6
Immunogen	Human CHTF18
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	CTF18, chromosome transmission fidelity factor 18 homolog (S. cerevisiae);CHTF18;C16orf41;C321D2.2;C321D2.3;C321D2.4;CHL12;Ctf18;RU VBL ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHTF18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.