



CLEC18A Antibody

Product Code	CSB-PA005521GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5D8T8
Immunogen	Human CLEC18A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	C-type lectin domain family 18, member A;CLEC18A;FLJ77955;MGC164804;MRCL;MRCL1 ;MRLP2
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLEC18A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.