



CLU Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA005595GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P10909 |
| Immunogen | Human CLU |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB,IHC |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Alias | clusterin;CLU;AAG4;APOJ;CLI;KUB1;MGC24903;SGP-2;SGP2;SP-40;TRPM-2; TRPM2 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | CLU |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |