



CNOT1 Antibody

Product Code	CSB-PA005663GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5YKK6
Immunogen	Human CNOT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	CCR4-NOT transcription complex, subunit 1;CNOT1;AD-005;CDC39;DKFZp686E0722;DKFZp686O168;FLJ36492;FLJ90644;KIAA1007;NOT1;NOT1H ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNOT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.