



COL6A2 Antibody

Product Code	CSB-PA005752GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12110
Immunogen	Human COL6A2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	Antigen Affinity purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COL6A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.