



COLEC11 Antibody

Product Code	CSB-PA005761GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWP8
Immunogen	Human COLEC11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Zebrafish
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	collectin sub-family member 11;COLEC11;CL-K1-I;CL-K1-II;CL-K1-IIa;CL-K1-IIb;CLK1;DKFZp686N1868;MGC3279 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	COLEC11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.