



# Recombinant Chlorocebus aethiops Serine/threonine-protein kinase 4 (STK4)

<b>Product Code</b>	CSB-MP022855DSU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A4K2Y1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	METVQLRNPP RRQLKKLDED SLTKQPEEVF DVLEKLGEGS YGSVYKAIHK ETGQIVAIKQ VPVESDLQEI IKEISIMQQC DSPHVVKYYG SYFKNTDLWI VMEYCGAGSV SDIIRLRNKT LTEDDIATIL QSTLKGLEYL HFMRKIHDI KAGNILLNTE GHAKLADFGV AGQLTDTMAK RNTVIGTPFW MAPEVIQEI YNCVADIWSL GITAIEMAEG KPPYADIHPM RAIFMIPTNP PPTFRKPELW SDNFTDFVKQ CLVKSPEQRA TATQLLQHPF VKSAKGVSI RDLINAMDV KLKRQESQQR EVDQDDEENS EEDMDSGTM VRAVGDEMGT VRVASTMTDG ANTMIEHDDT LPSQLGTMVI NTEDEEEEGT MKRRDETMQP AKPSFLEYFE QKEKENQINS FGKSVPGPLK NSSDWKIPQD GDYEFLLKSWT VEDLQKRLLA LDPMMEQEIE EIRQKYQSKR QPILDAIEAK KRRQQNF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	STK4
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase 4 EC= 2.7.11.1 Cleaved into the following 2 chains: 1. Serine/threonine-protein kinase 4 37kDa subunit Short name= 2. MST1/N 3. Serine/threonine-protein kinase 4 18kDa subunit
<b>Expression Region</b>	1-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.