



CPLX4 Antibody

Product Code	CSB-PA005896GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7G2
Immunogen	Human CPLX4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	complexin 4;CPLX4;CPX-IV;CPXIV;DKFZp686A0185;DKFZp686O0683;FLJ41190;MGC125769;MGC125783;MGC125784 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CPLX4
Target Details	This gene likely encodes a member of the complexin family. The encoded protein may be involved in synaptic vesicle exocytosis.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.