



CREB3L4 Antibody

Product Code	CSB-PA005952GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEY5
Immunogen	Human CREB3L4
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	cAMP responsive element binding protein 3-like 4; CREB3L4; AIBZIP; ATCE1; CREB3; CREB4; JAL; hJAL ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CREB3L4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.