



CRLS1 Antibody

Product Code	CSB-PA005985GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJA2
Immunogen	Human CRLS1
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	cardiolipin synthase 1; CRLS1; C20orf155; CLS1; GCD10; dJ967N21.6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CRLS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.