



# CRTC1 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA006000GA01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                           |
| <b>Uniprot No.</b>         | Q6UUV9  |
| <b>Immunogen</b>           | Human CRTC1   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human,Mouse,Rat   |
| <b>Tested Applications</b> | ELISA,WB  |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.    |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Alias</b>               | CREB regulated transcription coactivator 1;CRTC1;FLJ14027;KIAA0616;MECT1;TORC1;WAMTP1 ; |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | CRTC1   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.             |