



Recombinant Danio rerio Mitochondrial inner membrane protease ATP23 homolog (zgc:162885)

Product Code	CSB-MP026235DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4IGF3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAQSGAKAAD LSREPPGEQK PSPSSRQNEE DLGYNLYPER GSRLKKSTVE ESLFTFKHKC QLMLQFAMDT SPYAKLLLLGA MKSSGCTVFK DRHFSCEDCD GTVSGGFDA A TSQIVLCQNN IHQQAHMNRV VTHELIHAFD HCRAQVDWFS NYRHLACSEI RAANLSGDGS FINEFSRFNF GLRKHHQECV RGRALRSILA VRRVSREEAE RVVDEVFDSC FNDHAPFGRI PHSKKDAKFA HRDFENRDRY YANL
Source	Mammalian cell
Target Names	atp23
Protein Names	Recommended name: Mitochondrial inner membrane protease ATP23 homolog EC= 3.4.24.-
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.