



CTAG1A Antibody

Product Code	CSB-PA006125GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78358
Immunogen	Human CTAG1A
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	cancer/testis antigen 1A;CTAG1A;ESO1;LAGE2A ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	CTAG1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.