



CTNND1 Antibody

Product Code	CSB-PA006172GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60716
Immunogen	Human CTNND1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	catenin (cadherin-associated protein), delta 1;CTNND1;CAS;CTNND;KIAA0384;P120CAS;P120CTN;p120 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTNND1
Target Details	This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and althervative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined.