



CTR9 Antibody

Product Code	CSB-PA006178GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PD62
Immunogen	Human CTR9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Ctr9, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae);CTR9;KIAA0155;SH2BP1;TSBP;p150;p150TSP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTR9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.