



CTTNBP2 Antibody

Product Code	CSB-PA006208GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WZ74
Immunogen	Human CTTNBP2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	cortactin binding protein 2;CTTNBP2;C7orf8;CORTBP2;FLJ34229;KIAA1758;MGC104579;Orf4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTTNBP2
Target Details	This gene encodes a protein with six ankyrin repeats and several proline-rich regions. A similar gene in rat interacts with a central regulator of the actin cytoskeleton.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.