



# CUX1 Antibody

<b>Product Code</b>	CSB-PA006228GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q13948
<b>Immunogen</b>	Human CUX1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Alias</b>	cut-like homeobox 1;CUX1;CASP;CDP;CDP/Cut;CDP1;COY1;CUTL1;CUX;Clox;Cux/CDP;GOLIM 6;Nbla10317;p100;p110;p200;p75 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CUX1
<b>Target Details</b>	This protein is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.