



CX3CR1 Antibody

Product Code	CSB-PA006236GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49238
Immunogen	Human CX3CR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	chemokine (C-X3-C motif) receptor 1;CX3CR1;CCRL1;CMKBRL1;CMKDR1;GPR13;GPRV28;V28 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CX3CR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.