



# CXCL1 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA006239GA01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | P09341  |
| <b>Immunogen</b>           | Human CXCL1   |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Tested Applications</b> | ELISA,WB,IHC  |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.  |
| <b>Purification Method</b> | Antigen Affinity Purified   |
| <b>Isotype</b>             | IgG   |
| <b>Alias</b>               | chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha);CXCL1;FSP;GRO1;GRO $\alpha$ ;MGSA;MGSA-a;NAP-3;SCYB1 ; |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | CXCL1   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |