



CYB5B Antibody

Product Code	CSB-PA006310GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43169
Immunogen	Human CYB5B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	cytochrome b5 type B (outer mitochondrial membrane);CYB5B;CYB5-M;CYPB5M;DKFZp686M0619;OMB5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYB5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.