



Recombinant Aotus nancymae Serine/threonine-protein kinase 4 (STK4)

| | |
|--------------------------|---|
| Product Code | CSB-YP022855AUV |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A4K2W5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Aotus nancymae (Ma's night monkey) |
| Purity | >85% (SDS-PAGE) |
| Sequence | METVQLRNPP RRQLKKLDED SLTKQPEEVF DVLEKLGEGS YGSVYKAIHK ETGQIVAIKQ VPVESDLQEI IKEISIMQQC DSHHVVKYYG SYFKNTDLWI VMEYCGAGSV SDIIRLRNKT LTEDDIATIL QSTLKGLEYL HFMRKIHARDI KAGNILLNTE GHAKLADFGV AGQLTDTMAK RNTVIGTPFW MAPEVIQEIG YNCVADIWSL GITAIEMAEG KPPYADIHPM RAIFMIPTNP PPTFRKPELW SDNFTDFVKQ CLVKSPEQRA TATQLLQHPF VKSAKGVSIL RDLINAMDV KLKRQESQQR EVDQDDEENS EEDMDSGTM VRAVGDEMGT VRVASTMTDG ANTMIEHDDT LPSQLGTMVI NAEDEEEEGT MKRRDETMQP AKPSFLEYFE QKEKENQINS FGKSIPGLQ NSSDWKVPQD GDYEFLLKSWT VEDLQKRLLA LDPMMEQEIE EIRQKYQSKR QPILDAIEAK KRRQQNF |
| Source | Yeast |
| Target Names | STK4 |
| Protein Names | Recommended name: Serine/threonine-protein kinase 4 EC= 2.7.11.1 Cleaved into the following 2 chains: 1. Serine/threonine-protein kinase 4 37kDa subunit Short name= 2. MST1/N 3. Serine/threonine-protein kinase 4 18kDa subunit |
| Expression Region | 1-487 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |