



DCP1B Antibody

Product Code	CSB-PA006556GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZD4
Immunogen	Human DCP1B
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DCP1 decapping enzyme homolog B (S. cerevisiae);DCP1B;DCP1;hDcp1b ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DCP1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.