



DCUN1D5 Antibody

Product Code	CSB-PA006575GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTE7
Immunogen	Human DCUN1D5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	DCN1, defective in cullin neddylation 1, domain containing 5 (S. cerevisiae);DCUN1D5;FLJ32431;MGC2714 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DCUN1D5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.