



DCXR Antibody

Product Code	CSB-PA006577GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z4W1
Immunogen	Human DCXR
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dicarbonyl/L-xylulose reductase;DCXR;DCR;HCR2;HCR11;KIDCR;P34H;SDR20C1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DCXR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.