



Recombinant *Xenopus tropicalis* ELAV-like protein 4 (elavl4)

Product Code	CSB-MP007595XBF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4QNI8
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEWNGLKMI STMEPQVSNG PTSNTSNGPS SNSRNCPSM QTGAATDDSK TNLIVNYLPQ NMTQEEFRSL FGSIGEIESC KLV RDKITGT QFEEHFKDLA TGTKWKPLTE EGPIFGKGQS LGYGFVNYID PKDAEKAIN LNLRLQTKT IKVSYARPSS ASIRDANLYV SGLPKTMTQK ELEQLFSQYG RIITSRLVD QVTGVS RGVG FIRFDKRIEA EEAIKGLNGQ KPSGAAEPIT VKFANNPSQK TSQALLSQLY QSPNRRYPGP LHHQAQRFRL DNLLNMAYGV KRFSPITIDG MTSLVGMNIP GHTGTGWCIF VYNLSPDSDE SVLWQLFGPF GAVNNVKVIR DFNTNKCKGF GFVTMTNYDE AAMAIASLNG YRLGDRVLQV SFKTNKTHKS
Source	Mammalian cell
Target Names	elavl4
Protein Names	Recommended name: ELAV-like protein 4 Alternative name(s): Protein ElrD
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.