



Recombinant Danio rerio Phytanoyl-CoA hydroxylase-interacting protein-like (phyhipl)

Product Code	CSB-MP017949DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4QNW7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEVPRVGHSI SSPTSPCEDM MKNLSLDAIQ LCEREGNKSL DSGIAEMEEL PVPQNIKISN ITCDSFKICW DMDRSRKERI THYFIDLNKK ENKNANKFKH KDVPTKLVAK AVPLPMTVRG HWFLSPRTEY TVAVQTASKQ GDGDYAVSEW SEINEFCTAD YSTVHLTQLM EKAEVIAGRM LRFSVFYRNQ HKEYFDHARE AQNNKMLPSV KDNSGSHGSP ISGKLEGIFF SCNTEFNTGK PPQDSPYGRY RYQLPAEALF SPKTNLYFGD FYCMYTAYHY VILVIAPQGS PGDDFCKQRL PALDIANNRF MTCSEDEDGQ LVFHHAQDLI LEVIYDTPVD LSLGTVAEIS GHQLMSLSTV NAKKDPCKT CNISVGR
Source	Mammalian cell
Target Names	phyhipl
Protein Names	Recommended name: Phytanoyl-CoA hydroxylase-interacting protein-like
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.