



DDX43 Antibody

Product Code	CSB-PA006626GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXZ2
Immunogen	Human DDX43
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	DEAD (Asp-Glu-Ala-Asp) box polypeptide 43;DDX43;CT13;DKFZp434H2114;HAGE ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DDX43
Target Details	This protein is an ATP-dependent RNA helicase in the DEAD-box family and displays tumor-specific expression.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.