



DHRS4 Antibody

Product Code	CSB-PA006861GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTZ2
Immunogen	Human DHRS4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dehydrogenase/reductase (SDR family) member 4;DHRS4;DHRS4L2;FLJ11008;NRDR;SCAD-SRL;SDR-SRL;SDR25C1;humNRDR ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DHRS4