



DHX15 Antibody

Product Code	CSB-PA006870GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43143
Immunogen	Human DHX15
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DEAH (Asp-Glu-Ala-His) box polypeptide 15;DHX15;DBP1;DDX15;HRH2;PRP43;PRPF43;PrPp43p ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DHX15
Target Details	This protein is a putative ATP-dependent RNA helicase implicated in pre-mRNA splicing.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.