



DHX58 Antibody

Product Code	CSB-PA006884GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C10
Immunogen	Human DHX58
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DEXH (Asp-Glu-X-His) box polypeptide 58;DHX58;D11LGP2;D11lgp2e;LGP2 ;
Immunogen Species	Homo sapiens (Human)
Target Names	DHX58
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.