



DICER1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA006893GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UPY3 |
| Immunogen | Human DICER1 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | dicer 1, ribonuclease type III;DICER1;DCR1;Dicer;HERNA;KIAA0928 ; |
| Product Type | Purified Rabbit Anti Human,Mouse PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | DICER1 |
| Target Details | This gene encodes a protein possessing an RNA helicase motif containing a DEXH box in its amino terminus and an RNA motif in the carboxy terminus. The encoded protein functions as a ribonuclease and is required by the RNA interference and small temporal RNA (stRNA) pathways to produce the active small RNA component that represses gene expression. Two transcript variants encoding the same protein have been identified for this gene. |