



DIRAS1 Antibody

Product Code	CSB-PA006905GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95057
Immunogen	Human DIRAS1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DIRAS family, GTP-binding RAS-like 1;DIRAS1;Di-Ras1;FLJ42681;GBTS1;RIG ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DIRAS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.