



DMWD Antibody

Product Code	CSB-PA006979GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09019
Immunogen	Human DMWD
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dystrophia myotonica, WD repeat containing;DMWD;D19S593E;DMR-N9;DMRN9;gene59 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DMWD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.