



Recombinant Human Glycogen [starch] synthase, liver (GYS2)

Product Code	CSB-YP010079HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P54840
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLRGRSLSVT SLGGLPQWEV EELPVEELL L FEVAWEVTNK VGGIYTVIQT KAKTTADEWG ENYFLIGPYF EHNMKTQVEQ CEPVNDVAVRR AVDAMNKHGC QVHFGRWLIE GSPYVVLFDI GYSAWNLDRW KGDLEWACSV GIPYHDREAN DMLIFGSLTA WFLKEVTDHA DGKYVVAQFH EWQAGIGLIL SRARKLPIAT IFTTHATLLG RYLCAANIDF YNHLDKFNID KEAGERQIYH RYCMERASVH CAHVFTTVSE ITAIEAEHML KRKPDVVTN GLNVKKFSV HEFQNLHAMY KARIQDFVRG HFYGHLDLFDL EKTLFLFIAG RYEFNSKGD IFLESLRLN FLLRMHKS DI TVMVFFIMPA KTNNFNVELT KGQAVRKQLW DVAHSVKEKF GKKLYDALLR GEIPDLNDIL DRDDLTIMKR AIFSTQRQSL PPVTTHNMID DSTDPI LSTI RRIGLFNRT DRVKVILHPE FLSSTSPLL MDYEEFVRGC HLGVFPSY YE PWGYTPAECT VMGIPSVTTN LSGFGCFMQE HVADPTAYGI YIVDRRFRSP DDSCNQLTKF LYGFCQSR R QRRIQRNRTE RLSDLLDWRY LGRYYQHARH LTLSRAFPDK FHVELTSPPT TEGFKYPRPS SVPPSPSGSQ ASSPQSSDVE DEVEDERYDE EEEAERDLN IKSPFSLSHV PHGKKKLHGE YKN
Source	Yeast
Target Names	GYS2
Protein Names	Recommended name: Glycogen [starch] synthase, liver EC= 2.4.1.11
Expression Region	1-703
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein, liver glycogen synthase, catalyzes the rate-limiting step in the synthesis of glycogen - the transfer of a glucose molecule from UDP-glucose to a terminal branch of the glycogen molecule. Mutations in this gene cause glycogen storage disease type 0 (GSD-0) - a rare type of early childhood fasting hypoglycemia with decreased liver glycogen content.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.