



DNAJB1 Antibody

Product Code	CSB-PA007003GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25685
Immunogen	Human DNAJB1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DnaJ (Hsp40) homolog, subfamily B, member 1;DNAJB1;HSPF1;Hdj1;Hsp40;Sis1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.