



DNAJB6 Antibody

Product Code	CSB-PA007012GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75190
Immunogen	Human DNAJB6
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DnaJ (Hsp40) homolog, subfamily B, member 6;DNAJB6;DJ4;DKFZp566D0824;DnaJ;FLJ42837;HHDJ1;HSJ-2;HSJ2;MGC1152;MGC117297;MRJ;MSJ-1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJB6
Target Details	This gene encodes a member of the DNAJ protein family. DNAJ family members are characterized by a highly conserved amino acid stretch called the J-domain and function as one of the two major classes of molecular chaperones involved in a wide range of cellular events, such as protein folding and oligomeric protein complex assembly. This family member may also play a role in polyglutamine aggregation in specific neurons. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been fully described.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.