



# DNAJC16 Antibody

|                            |                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA007023GA01HU                                                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                             |
| <b>Uniprot No.</b>         | Q9Y2G8                                                                                                    |
| <b>Immunogen</b>           | Human DNAJC16                                                                                             |
| <b>Raised In</b>           | Rabbit                                                                                                    |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                           |
| <b>Tested Applications</b> | ELISA,WB                                                                                                  |
| <b>Storage Buffer</b>      | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.                     |
| <b>Purification Method</b> | Antigen Affinity purified                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                       |
| <b>Alias</b>               | DnaJ (Hsp40) homolog, subfamily C, member 16;DNAJC16;DKFZp686G1298;DKFZp686N0387;DKFZp781I1547;KIAA0962 ; |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                                            |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                      |
| <b>Target Names</b>        | DNAJC16                                                                                                   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                               |