



DNAJC4 Antibody

Product Code	CSB-PA007036GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NNZ3
Immunogen	Human DNAJC4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	DnaJ (Hsp40) homolog, subfamily C, member 4;DNAJC4;DANJC4;HSPF2;MCG18;MGC19482;MGC57189;MGC71863 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DNAJC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.