



DOCK3 Antibody

Product Code	CSB-PA007099GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZD9
Immunogen	Human DOCK3
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	dedicator of cytokinesis 3; DOCK3; KIAA0299; MOCA; PBP ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DOCK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.