



DPP3 Antibody

Product Code	CSB-PA007138GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY33
Immunogen	Human DPP3
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA
Form	Liquid
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DPP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.