



DTX3L Antibody

Product Code	CSB-PA007224GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDB6
Immunogen	Human DTX3L
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	deltex 3-like (Drosophila);DTX3L;BBAP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DTX3L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.