



# Recombinant *Magnaporthe oryzae* Flap endonuclease 1 (FEN1)

<b>Product Code</b>	CSB-EP008585MPM
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A4QS18
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Magnaporthe oryzae</i> (strain 70-15 / ATCC MYA-4617 / FGSC 8958) (Rice blast fungus) ( <i>Pyricularia oryzae</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MGIKQLFTII KENAPAAIKT GEIKNQFGRK VAIDASMSIY SFLIAVRSNG            EMLTNEDGQT TSHLMGMFYR TLRMVDNGIK PLYVFDGAPP KLKSGELARR            YQRKQEALLEG LEEARETGTA EDVEKFSRRT VRVTREHNEE CRQLLKMGI            PYIIAPTEAE AQCAVLARAG KVFAAASEDM DTLCFDSPIL LRHLTFSEAR            KEPIQEIHVD KVLEGLDMDR KQFVDLCILL GCDYLDPIPK VGPSTALKLI            REHGSLETIV EKMKKGELKY TVPEDWPFED ARDLFFNPAV HPADHPDCNF            KWEKPDVEGL VKYLVTEKGF SEDRVRPGAL RLEKALGTSQ QQRLEGFFKP            VARTAEQKA HKRKLEVKAE EAKKKLKA EKAKAKARP RGTA         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	FEN1
<b>Protein Names</b>	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
<b>Expression Region</b>	1-394
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.