



EEF2 Antibody

Product Code	CSB-PA007434GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13639
Immunogen	Human EEF2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	eukaryotic translation elongation factor 2;EEF2;EEF-2;EF2 ;
Product Type	Purified Rabbit Anti Homo sapiens PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EEF2
Target Details	This gene encodes a member of the GTP-binding translation elongation factor family. This protein is an essential factor for protein synthesis. It promotes the GTP-dependent translocation of the nascent protein chain from the A-site to the P-site of the ribosome. This protein is completely inactivated by EF-2 kinase phosphorylation.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.